

USA DENTAL SUPPLIES INC

Product Price List (US) — Effective March 2026

Delivering Quality and Price




Authorized by Li Liu, President

EMAIL
usds2030@gmail.com

PHONE
(470) 846-1775

WEBSITE
usadental.shop

Picture	Product Name	Manufacturer	Description	Retail Price (USD)	March 2026 Promotion
	Flight A6 Operatory Package	Flight Dental Systems	Complete dental operatory package — hybrid hydraulic-electromechanical patient chair with Torch LED light, delivery system, cuspidor, and full assistant instrumentation (HVE, SE, air/water syringe). 400 lb capacity, 12 upholstery colors, North American or Continental delivery style.	\$14,664.00	\$11,731.00
	Dental Camera Recorder System (DCR-1)	Medfruition	Full HD 1080P dental camera with 30x optical zoom, Sony sensor, 8,000–30,000 lux LED illumination (Rae90). Supports live broadcasting, remote consultations, and chairside education. Records MP4; plug-and-play setup.	\$9,600.00	\$7,600.00
	Runyes 3DS V5 Intraoral Scanner	Runyes Medical Instrument Co.	Lightweight (160 g) intraoral scanner with <10 µ single-crown accuracy, 25 mm depth of field, motion-sensing gesture control. Open STL/PLY/OBJ output — compatible with Exocad, 3Shape and all major CAD/CAM platforms.	\$7,400.00	\$6,000.00
	iRay Pluto Intraoral Sensor (Size 1.5 & 2)	iRay	20 µ pixel pitch CMOS sensor with direct CsI scintillator. High resolution, low dose digital radiography. AED trigger, fast and easy workflow.	\$2,888.00	\$2,500.00
	iRay Pluto Intraoral Sensor (Size 0 & 1)	iRay	20 µ pixel pitch CMOS sensor with direct CsI scintillator. High resolution, low dose digital radiography. AED trigger, fast and easy workflow.	\$2,665.00	\$2,300.00
	Dental Burs — Carbide & Diamond	Yepdent	Complete range of carbide and diamond burs in FG and RA shanks. 150+ carbide variants, 180+ diamond variants across all shapes and grits, including a Zirconia Diamond line. Packed in boxes of 10. See separate bur price sheet for full SKU list.	from \$19.00 / box	from \$15.00 / box
	Biosis Healing Membrane	Biosis Healing	Resorbable collagen barrier membrane for guided bone regeneration (GBR) and guided tissue regeneration (GTR). Resorbs naturally — no second surgery needed. Compatible with all bone graft materials. Five sizes: 15x15 mm to 30x40 mm.	from \$69.00	from \$41.00

	<p>Surgical Drive System (ES70)</p>	<p>Jinme Medical</p>	<p>All-in-one intelligent surgical drive system for implant placement, tooth extraction, restoration, and oral surgery. Includes 4 hand-pieces (1:4.2, 1:5, 1:1, 40:1), peristaltic pump, 0–200,000 RPM. No fixed water or air supply required.</p>	<p>\$4,600.00</p>	<p>\$4,600.00</p>
	<p>Endo Motor</p>	<p>Bondent / Bomedent</p>	<p>Cordless endodontic motors — 4 models: Endo Wise, Endo Wise Plus (with apex locator), iRoot Pro, and Marc III. Multiple reciprocating modes, large OLED display, accurate torque control, extended battery life.</p>	<p>\$980.00</p>	<p>\$980.00</p>
	<p>Dental Operating Microscope</p>	<p>Bondent</p>	<p>Two models: DOM 4000D 4K (G4) with built-in 16MP 4K camera, Wi-Fi streaming, F=180–460mm optics, 1,650mm arm. G4 Pro with Sony 4K sensor, 12MP stills, apochromatic German optics, orange/green/fluorescence filters, electronic focus, and Bluetooth pedal control.</p>	<p>Contact for pricing</p>	<p>Contact for pricing</p>

Ready to Order or Need a Quote?

EMAIL
usds2030@gmail.com

PHONE
(470) 846-1775

WEBSITE
usadental.shop